



The University of Sheffield

Risk Assessment Form

PERSONS AT RISK : () Employees () Contractors () Public () Visitors () Others			Reference No:
Risk: (H) High (M) Medium (L) Low (O) No Risk.		Environment: Ground level laboratories with good access / egress	
TASK or ACTIVITY: Operation and maintenance of 500 MHz & 600 MHz NMR Spectrometers		INITIAL RISK RATING	EXISTING CONTROL/PROPOSED CONTROL MEASURES
SIGNIFICANT HAZARD	RISK		
Accidental release of significant volumes of liquid or gaseous cryogen (Helium or Nitrogen)	Asphyxiation	M	1. Low level oxygen deficiency sensing & alarm system] 2. Emergency ventilation triggered by low oxygen levels - 100% fresh air] 3. Restricted access arrangements - authorised staff only.] 4. All oxygen deficiency sensors checked weekly, calibrated 6 monthly and new sensor installed every 12-18 months.] L 5. Emergency procedure for action in event of gas alarm implemented.]
Contact with liquid cryogen during refilling sequence	Cold burns to exposed skin	M	1. All operators must wear laboratory coat, suitable, dry leather gloves, eye protection, suitable leg protection and stout shoes.] 2. Refilling operations accord with Safe Operating Procedure] L 3. All operators to be trained in all relevant Safe Operating Procedures]
Strong magnetic field	Injury from moving ferromagnetic objects	M	1. All operators and visitors required to remove all ferromagnetic objects prior to entry to laboratory.]
	Pacemaker failure		2. Warning signs posted on access door to laboratory] L
Comments:			Overall Risk:
Oxygen deficiency sensing equipment checked weekly by introducing oxygen deficient source to test detection and warning systems.			L
Additional References, Tasks etc			
Safe Operating Procedures for Refilling NMR unit with Liquid Helium and Liquid Nitrogen			
All users of NMR's to be familiar with the document "General Safety in the Use of Cryogenic Liquids" on the Safety Services web site.			
Undertaken By:			
Other Persons Consulted:			
Date:			Revision Date: